

OPTIMAL PERFORMANCE MADE SIMPLE



SS-501SX

Microprocessed Electrosurgical Unit

The SS-501SX Microprocessed Electrosurgical Unit adds value to the OR by combining the basics of an electrosurgical generator with specific safety and control features.


Medtronic
Further, Together

CONTROL

Achieve desired clinical performance by choosing between the twenty monopolar and bipolar modes available.

SAFETY

Built-in safety features, including patient return electrode monitoring, to block energy delivery in out-of-specification impedances between the pad and the patient.



PERFORMANCE

EZcut system provides cutting performance for monopolar and bipolar cut modes, while the initialization auto-testing ensures designed performance in everyday use.

SIMPLE

Remote function simply enable power out put regulation from hand switching accessories in the surgical field.



SS-501SX BASIC TECHNICAL SPECIFICATIONS

Product Dimension (D X W X H)
397 x 300 x 163 mm

Product Weight
5.45 kg

Rated Supply Voltage Range
100 to 240 VAC with auto voltage selection

Rated Frequency Range
50 / 60 Hz

Max. Rated Input
1 kVA

Liquid Ingress/Spillage — SS-501SX
IP21

Liquid Ingress / Spillage — Footswitches
IPX7

Max. Power Output
400 W

Power Adjustment Range Increment
1 Watt (and Microbipolar 0.5 Watt)

Power Output Tolerance
15% or 5 Watts (whichever is higher)

Basic Frequency of the Unit
390 kHz sinusoid

| Ordering Code | Trade Model | Product Description |
|---------------|-------------|-------------------------|
| 06725 | SS-501SX | Argon Plasma Coagulator |

ACCESSORIES

| | | |
|----------------------|-------|---|
| 06523 | FS-16 | Footswitch, double pedal, Monopolar 1, for SS-501SX |
| 06524 | FS-17 | Footswitch, double pedal, Monopolar 2, for SS-501SX |
| 09922 | FS-25 | Footswitch, single pedal, Bipolar, for SS-501SX |
| 09921 | FS-24 | Footswitch, double pedal, Bipolar, for SS-501SX |
| Check options | - | Power cords |

For more information, contact your Medtronic energy representative

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 07/2021 – US-SE-2100117 – [WF # 4930843]

550 Marechal Mascarenhas de Moraes
14095-120, São Paulo, Brazil

rs.WEM@medtronic.com

www.wem.com.br

+55 (16) 3512-4600

Medtronic